


PI-Max3


**Global handles and fast data collection
examples using the PI-Max3**

Open Toolkit and Camera

Open ExOpenGlobalCam.vi which can be found in your :
SIToolKit\LabView\Examples\xx_ver\Fast Examples\Open-Close Hardware
folder. This only needs to be run once when you enter LabView. Note camera
should be turned on and attached to the computer before this is run.




If the camera opens successfully you will get a **green** light in the LED on the front panel. There is now a global variable that holds the camera handle and camera open status for other programs to use.




Open Pulser

Run the ExOpenGlobalPulser.vi

This opens the Pulser and stores the pulser handle into a global variable to be used by other VI's




. If the pulser is opened successfully the LED will turn **Green**




Setup PI-Max3 For Internal triggering.


This is good to try for a sanity check to make sure the system is working correctly as far as software and hardware are concerned since it is self contained and does not rely on outside triggering.



Run ExDataCollectPulserFocus.vi –
this will verify that the hardware is working. Wait until the
ready to start led turns green then hit the start button.




The focus program should display an image that is updated. Also on the back of the of the PI-Max3 the green gate indicator light should flash on and off.



High speed data collection:


Run the ExDataCollectionPulserNFrame.vi This program is given a number of images to collect. It allocates in RAM enough buffer space (below Labview in C++) to hold this amount and then uses only a 4 byte handle in LabView for the actual data collection loop.

First wait for the ready to start LED is green then hit the start button. The data will be stored without display. When it is done a beep will sound and the program will then store the data into the file nframedata.




PI-Max3 in external triggering experiments.

Run ExSetupExtSuperSyncro.vi (in place of running ExSetupIntSuperSyncro.vi) to setup external triggering with the PI-Max3. After this both the ExDataCollectPulserFocus.vi and the ExDataCollectPulserNFrame.vi can be run to collect data.



Before leaving LabView the software should be unloaded to deallocate buffers and close down the hardware. With a PI-Max system, first close down the Pulser section by running ExCloseGlobalPulser.vi



Then close out the camera, which will also close out the toolkit

